


Overview

Reliable, Reproducible, and Research-Proven

AstroPure Bead is a bead-based nucleic acid purification solution designed to deliver consistent performance across a wide range of genomic applications. From sequencing to PCR, AstroPure helps maximize DNA recovery and streamline workflows, ensuring clean, high-quality output at every step.



Why choose AstroPure Bead?

Outstanding Recovery Optimised for amplicons above 100bp, AstroPure retains more of the target DNA.

Trusted Across the Industry Cited in over 15,000 studies and included in more than 200 commercial library prep kits from industry leaders such as Illumina, IDT, and Thermo Fisher.

Efficient Impurity Removal Effectively removes short fragments, primers, salts, and enzymes that may compromise downstream analysis.

Flexible Compatibility Suitable for both manual and automated workflows with identical protocols.

Consistent Results Lot-to-lot consistency and high-quality output support high-stringency analyses.

Applications

NGS Library Preparation Removes enzymes, unreacted components, and low-molecular-weight contaminants for optimal sequencing accuracy.

qPCR/ddPCR Workflows Provides clean DNA templates free from inhibitors, supporting sensitive quantification and detection.

Sanger Sequencing Improves chromatogram clarity by removing excess primers and dNTPs post-PCR.

Microarray Preparation Delivers high-purity DNA for precise hybridization and signal interpretation.

Post-PCR Purification Supports cloning, restriction enzyme digestion, and labelling through the removal of residual reagents.

Superior Performance in NGS Clean-Up

With an optimal bead-to-sample ratio (1.8x), AstroPure Bead provides an effective balance between performance and cost-efficiency. DNA yields consistently fall within the ideal range (80-95%) for next-generation sequencing, reducing data loss and the need for repeat experiments.

AstroPure Bead vs Other kits

Surpasses 10+ alternative products in DNA recovery rate comparisons.

- Maintains procedural simplicity while high analytical stringency.
- Minimises background interference in Sanger sequencing compared to column-based purification methods.

Compatible Sample Types

- Fresh/frozen tissue
- FFPE (formalin-fixed, paraffin-embedded) tissue
- Cultured cells
- Blood- and plasma-derived nucleic acids
- Low-input or degraded DNA

Compatible Workflows

- Manual bench-top processing
- High-throughput automation (e.g. 96-well format)
- Downstream applications such as cloning, ligation, and hybrid capture

